

ORION®

ORION ELECTRONICS LTD (EUROPE)  

SPA FOOT WARMER WITH MASSAGE USER MANUAL

MODEL: ORION OWP20F

You are now able to have full benefit of the massage that has proven to therapeutically regenerate tired, knotted muscles, relieve the tension of stressful living and provide an overall dimension of quality life.

IMPORTANT SAFETY INSTRUCTIONS

Read all instruction before using the cushion.

DANGER – To reduce risk of electric shock

- Always unplug from the electrical outlet immediately after using and before cleaning.
- Do not reach for an appliance that has fallen into water. Unplug immediately.
- Do not use while bathing or in a shower.
- Do not place or store appliance where it can fall or be pulled into a tub sink. Do not place in or drop into water or other liquid.

WARNING – To reduce risk of burns, fire, electric shock, or injury to persons:

- An appliance should never be left unattended when plugged in. unplug from outlet when not in use and before cleaning.
- To disconnect, turn control to OFF position and then remove plug from outlet.
- Close supervision is necessary when this appliance is used by, on or near children, invalids or disabled persons.
- Use this appliance only for its intended use as described in this manual. Do not use any attachment not recommended by the manufacturer.
- Never operate the appliance if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return this appliance to a service center for examination and repair.
- Do not carry this unit by the power cord.
- Do not operate under blanket. Excessive heating can occur and cause fire, electric shock or injury to persons.
- **Continued operation for more than 15 minutes is not recommended. Extensive use could lead to the produce's overheating. In the event of overheating, discontinue use and allow the unit to cool before operating.**
- Disconnect the power cord by holding the plug attached to the cord, not by using the cord itself.
- Never drop or insert any object into an opening.
- Do not operate where aerosol (spray) product is being used or where oxygen is

being administered.

- Do not operate this massager near an electrical oven or any heat source.

CAUTIONS

- Do not exceed 15 minutes in any one area of the body.
- Never use this massager as a substitute for medical attention. Persistent ailments and symptoms could indicate a more serious condition.
- Never use this massager directly on open wounds or bruised areas.
- Not to be used by children below 12 years of age.
- Do not use by pregnant, high blood pressure and heart attack.
- Never use while sleeping.
- No not use on or near eyes or other highly sensitive areas.

OPERATING INSTRUCTIONS

- Start to massage by move the buttons on controller.
- To turn on motor, move the MASSAGE button to HI or LO position. To turn off, move the MASSAGE button to OFF position.
- To turn on heat, move the HEAT button to ON position. To turn off, move the HEAT button to OFF position.



MAINTENANCE AND CLEANING

- The unit is maintenance free. If a fault does occur, just check whether the mains power adapter is properly connected.
- If the unit becomes faulty, do not try to repair it yourself.

TECHNICAL DATA

Model: ORION OWP20F
Size: 28 x 30,5 x 24cm
Material: polyester 100%
Net weight: 0,56kg
Colour: black

Cord length of controller: 75cm
Cord length of adapter: 180cm
AC 230V~50Hz
DC 12V, 600mA, 7.2VA
Protection class: II